



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Water Baths Thermostatic

Product Code :
NMCBS11063



Description :

Water Baths Thermostatic We are offering a wide range of water baths thermostatic that are fabricated from top quality materials as per the international quality standards. These water baths are known for their high efficiency, robust construction and corrosion resistant. Some of the features are: •Double walled construction •Inner made of stainless steel 304 & Outer made of mild steel with powder coating finish. •Temperature controlled by Thermostat. •Temperature Range : 5°C above ambient to 95 °C •6 or 12 holes 7.5cm dia with set of concentric rings for each hole. Water Baths Thermostatic, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,