



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Autoclaves (Steriliser) Portable

**Product Code :**  
NMCBS11067



**Description :**

Autoclaves (Steriliser) Portable We offer a wide range of autoclaves that are sterilize and are portable one. These autoclaves can work on gas or electricity. The container is made of aluminum or stainless steel. The Lid is fitted with pressure gauge & safety valve which is adjusted to blow off the steam at 12 psi maintaining constant pressure of 15 psi. Further, the technical specifications of these incubators are: Size MOC Cap (Liters) Size (cm) (Dia x Ht) AEI-14a Aluminium 20 30 x 30 AEI-14B St.Steel 20 30 x 30 Autoclaves (Steriliser) Portable, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,