

Product Name :
Muffle Furnace**Product Code :**
NMCBS11079**Description :**

Muffle Furnace The most cost effective bench-top high temperature Muffle Furnace made by our company is widely used in hospitals, clinics, research centres and other institutions. These muffles are rectangular in shape and light in weight. They are completely ceramic muffle with heat fundamentals of excellent kanthat wire. Few Features of Muffle furnaces are: Pyrometer Two indicating lamps Temperature Maxm 1000°C / 1200°C Thermocouple lead wire Fibre wool insulation 230V AC Mains 50HZ Working: 940°/1100° controlled by Digital temperature controller Some of the technical specifications are: Size Muffle Chamber Depth x Height x Width (cms) Rating Watts AEI-18a 22.5 x 10.0 x 10.0 1600 AEI-18b 25.0 x 12.5 x 12.5 2000 AEI-18c 30.0 x 12.5 x 12.5 3000 AEI-18d 30.0 x 15.0 x 15.0 3500 AEI-18e 45.0 x 22.5 x 22.5 6000 Muffle Furnace, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,