

**Product Name :**  
Incubator**Product Code :**  
NMCBS11081**Description :**

Incubator We offer incubator of memmert that are widely use in hospitals, clinics and other institutions. These incubators are known for their robust body, durability and long service life. Some of the features of our range of memmert incubator are: •Double walled construction •Exterior made of mild steel powder coated and interior made of stainless steel 304 (STD Model) or Exterior and interior made of stainless steel 304 (GMP Model). •80/20 Ni Heating elements placed at bottom and two vertical sides for uniform heating. •Arrangement for L shapes Thermometer facing top front side. •Fan arrangement for forced air circulation •Digital Temperature controller •Temperature Range : 5°C above ambient to 80°C •Built in glass window in the front door is provided for inner observation. •Power supply: 230 V AC Mains 50 HZ. Further, the technical specifications of these incubators are: Models Chamber Watts Shelves H x W x D (cms) AEI-18a 30 x 30 x 30 250 1 AEI-18b 35 x 35 x 35 300 2 AEI-18c 45 x 45 x 45 500 2 AEI-18d 60 x 45 x 45 600 2 AEI-18e 60 x 60 x 60 750 3 AEI-18f 90 x 60 x 45 900 3 AEI-18g 90 x 60 x 60 1000 3 Incubator, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers