



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Environmental Test Chamber

Product Code :
NMCBS11082



Description :

Environmental Test Chamber Double walled construction Exterior made of stainless steel 304 Interior made of stainless steel 316 Puff insulation (75 mm) Full length inner glass door for clear inner view Fan arrangement for forced air circulation through out the chamber for uniform temperature Pid temperature controller cum indicator Cooling system : hermetically sealed compressor unit coupled with evaporation coil & condenser. Temperature range : 10°C to 60°C Resolution: +/-1°C Humidity range : 40%rh to 95%rh Resolution: +/-3%rh Environmental Test Chamber, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,