



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Mortuary Chamber

Product Code :
NMCBS11083



Description :

Mortuary Chamber Technical Specification Mortuary Storage System is designed for storing cadaverous under cool condition to prevent decomposition. Double Walled with Mild Steel panels for outer surface and SS inner chamber with PUF insulation in between both the walls. Front opening, hinged insulated doors lined with magnetic gasket, handle and lock arrangements with keys in duplicate for individual dead bodies. Trays to carry dead bodies with handles at both ends, travel rollers fitted with ball bearing for easy movement. Exterior surface is chemically treated, antirust coated and duly furnished with hammer tone paint. CFC free Refrigeration unit is kept outside the chamber and fixed at rear top. Temperature of 2°C to 6 °C is maintained by hermetically sealed refrigeration unit and fitted with 'EMERSON' make compressor. Temperature is controlled by an automatic electronic solid state Digital Temperature controller cum Indicator. Supplied complete with trays, cord & plug. Suitable to work on 220 V Single Phase 50 HZ Ac supply. Mortuary Chamber, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,