



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Heating Ovens

Product Code :
NMCBS11088



Description :

Heating Ovens We are manufacturers of Heating Ovens. These are known for their robust construction, low maintenance and high efficiency. Some of the features are as follows: •Double Walled Construction •Exterior & Interior made of stainless steel 304/316 (GMP model) or Exterior made of Mild Steel with powder coating finish and interior made of stainless steel 304 (STD Model) •80/20 Ni Heating Elements placed at bottom •Temp. Range : 5°C above ambient to 300°C 2°C •230V AC 50HZ Mains Supply. •Digital Temperature Controller. •Fan arrangement for forced air circulation. Heating Ovens, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE