

Product Name :
Dehumidifiers**Product Code :**
NMCBS11098**Description :**

Dehumidifiers offered are designed to efficiently remove excess moisture from atmosphere of laboratory, work place or from storage facility. Some of the end application areas where these are used in include research laboratories, pharmaceutical industries, seed storage facility and other areas that require dry air or low humidity atmosphere for different applications. Here, the 1 ton dehumidifier being offered is designed for floor area of 100 sq ft and 1.5 ton dehumidifier is designed for floor area of approximately 150 sq ft. With the system working on 220 volts, single phase and 50Hz ac supply, these provide for ideal freezing environment for advanced medical biotechnological as well as industrial applications. Some of the standard features include optimum temperature of up to -35.0°C, micro-controller based control, bright red LED display and temperature stability with reduced energy consumption through PUF insulation. Dehumidifiers, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,