



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Soxhlet Extraction Heater

**Product Code :**  
NMCBS11099



**Description :**

Soxhlet Extraction Heater provided finds usage where several projects require to be run simultaneously or multiple runs of same process is needed. Here, the construction of the equipment comprises six 4 1/2" cast aluminum plates that having individual heat control. Further, these can have same temperatures for all six plates else temperatures can also be adjusted individually. For optimum performance, all hot plates have 300 W heaters with temperatures of equipment ranging from normal to 300°C. Here, the equipment also comes with pilot lamps for individual plates and lamp for power supply. Further, the surface damage can be prevented through use of rubber buffers. The other operational features of these extraction heaters comprise utility rack comprising of 12 clamps, 12 clamp holder and 4 rods with the system running on 220 to 230 AC volt option Soxhlet Extraction Heater, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,