



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Low Temperature Water Bath

**Product Code :**  
NMCBS11100



**Description :**

Low Temperature Water Bath offered comprise standard and customized Low Temperature water Baths that are developed for perfectly matching up to the challenging work process demands of scientists while they handle and specialized research projects. Here, the systems are widely used for day bacteriological as well as various laboratory applications requiring incubation, general tests and laboratory procedures like Wassermann and Kilmer tests. Further, the Low Temperature Water Baths also finds usage in different tissue culture applications, growth observation studies, enzyme reaction studies, fermentation analysis as well as other general and specialized laboratory applications. Some of the end usage areas where these Low Temperature Water Baths are supplied to include Research laboratories, Defence installations, Educational institutes as well as R & D labs of national and international repute. For superior performance, these Refrigerated Water Bath comprises double walled convection heated unit with outer wall made of thick PCRC steel sheets that is pre-treated with primer for rust proofing as well as painted with air drying paint / stove enamel / powder coating so as to deliver optimum service life standards.



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,