


Product Name :
Thermostatic Water Bath**Product Code :**
NMCBS11112**Description :**

Thermostatic Water Bath We offer our clients premium quality Thermostatic Water Bath. These are also known as serological water Bath which works on the temperature ranging from ambient temperature to 80° C. Thermostatic Water Bath temperature is controlled by the thermostatic with accuracy of $\pm 0.5^{\circ}\text{C}$. These products are double walled inside the stainless steel and outside are painted with mild steel in epoxy powder coating. Our products consists of two pilot lamps and temperature control knob and on or off switch. These works on 220/230V AC and are supplied with or without stirring arrangement without any racks and thermometer. Our product has a lid of water bath and is made of stainless steel sheet. Optional: · L-shaped thermometer (fitted with water bath) · Imported thermostat "EGO" or "JUMO" German with an accuracy of $\pm 0.5^{\circ}\text{C}$ · Digital temp. indicator-cum-controller · Microprocessor PID digital temp. indicator-cum-controller · Stirrer with 1/20 H.P. High speed motor with stainless steel stirring rod and blade Specification: Size inside chamber (Stainless Steel) Capacity (Approx.) 300x200x175 mm (suitable for 2 racks) 14 ltrs. 375x300x175 mm (suitable for 4 racks) 18 ltrs. 455x300x175 mm (suitable for 6 racks) 25 ltrs. 605x300x175 mm (suitable for 8 racks) 32 ltrs. Thermostatic Water Bath, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,