



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Refrigerated Liquid Bath

Product Code :
NMCBS11113



Description :

Refrigerated Liquid Bath Being a customer focused organization, we are engaged in offering a wide and exotic assortment of Refrigerated Liquid Bath. These products has a standard construction with the inner chamber manufactured of stainless steel sheet with digital temp. controller -cum-indicator. Refrigerated Liquid Bath temperature can be controlled ranging from minus (-) 20°C to plus (+) 20°C. Our products are made made as per the German standards and availed as per the needs and specifications of our clients. These products are highly qualitative as these are manufactured from the qualitative raw material procured from the trusted vendors. Features: Refrigerated Liquid Bath Sturdy construction Refrigerated Liquid Bath high performance Refrigerated Liquid Bath longer service life Specification: Chamber Size Capacity (Approx.) 450x300x300 mm 42 ltrs. Refrigerated Liquid Bath, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,