

**Product Name :**  
Blood Bank Refrigerator**Product Code :**  
NMCBS11116**Description :**

Blood Bank Refrigerator With our expertise in this domain, we are able to offer a superior quality range of Blood Bank Refrigerator. These products offered by our organization are the most versatile unit for the blood bank which is mounted on the sturdy angle frame. Blood Bank Refrigerator inner chamber is made up of stainless steel sheet and the outer surface is manufactured of the mild steel sheet duly powder coated paint. These products are provided for easily mobility. Our products are maintained in 4°C to 6°C with forced air circulation ensuring the uniform temperature in order to store the blood bottles or blood bags. Blood Bank Refrigerator temperature of these products are controlled by the microprocessor PID digital temperature controller-cum-indicator having a triggers the audio visual alarm during the fluctuation permitting the tolerance. These units also has door operated alarm and power supply failure alarm for security. Our products consists of efficient refrigeration system & hermetically sealed compressor for high performance maintaining the uniformity of temperature and also easy viewing chamber. Available in following capacities: · 85 liters (50 blood Bags) · 165 liters (110 blood Bags) · 200 liters (200 blood Bags) · 280 liters (250 blood Bags) · 400 liters (360 blood Bags) Optional: · Chart type temperature recorder · Voltage Stabilizer as recommended Blood Bank Refrigerator, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

