



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Vaccine Storage Cabinet

Product Code :
NMCBS11117



Description :

Vaccine Storage Cabinet Clients can avail from us a wide range of Vaccine Storage Cabinet. The inner chamber is made of stainless steel external mild steel sheet body of rigid steel frames. These are insulated thoroughly providing with face fitting doors which is complete with self-closing hinges lockable fasteners. Our products are fully sealed for avoiding the spilling and contamination from external environment. These works under the temperature ranging from 4°C to 6°C and are available in various sizes for keeping vaccines bottles. We have a drawer with clear perspex front and index card holders. Vaccine Storage Cabinet temperature is controlled and displayed by microprocessor PID digital temperature controlled-cum-indicator. Our products have a audio visual alarms which are triggered during the exceeding fluctuations permitted tolerance level and are suitable to work on 220V, single phase, 50Hz, AC supply. Available in following size: · 450x400x850 mm · 600x675x850 mm
Optional: · Chart type Temperature recorder · Voltage Stabilizer(Recommended) Vaccine Storage Cabinet, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,