

**Product Name :**  
Electrical Incinerator**Product Code :**  
NMCBS11118**Description :**

Electrical Incinerator Our clients can avail from us an exquisite range of Electrical Incinerator. These products offered by our organization are electrically or oil fired incinerator only first two sizes and can be supplied fully electrically operated with the other sizes shall be supplied oil fired (LDO). Our products are suitable for incineration of pathological waste. Electrical Incinerator product also includes cost of the chimney. These products are made as per the German standards. Electrical Incinerator Specification: · 5-8 Kg · 10-15 Kg · 25-35 Kg Electrical Incinerator, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)