

**Product Name :**  
Reynolds Apparatus**Product Code :**  
NMCBS11124**Description :**

Reynolds Apparatus Our organization is considered as one of the most trusted names of the industry engaged in offering genuine quality Reynolds Apparatus to the customers. These apparatuses are designed in compliance with the international quality standards. Moreover, our quality controllers check the entire range on certain standard parameters. Some Of The Quality Attributes Are As Follows: Material: Stainless steel Heat water tank capacity: 40 ltr Flow: Using measuring cylinder Reynolds Apparatus, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)