

**Product Name :**  
Mass Transfer Lab Equipment

**Product Code :**  
NMCBS11139



**Description :**

Mass Transfer Lab Equipment We are instrumental in offering superior quality Mass Transfer Lab Equipment to the esteemed customers, as per their specific demands. These mass transfer labs are designed in such a manner that our clients' exact requirements get fulfilled by us. In addition to this, we have gained expertise in offering customized products to the clients. Application Area Of Product: Solid-Liquid Extraction (Packed Bed Type) Absorption in Wetted Wall Column Packed Distillation Column Bubble Cap Distillation Column SIMPLE STEAM Distillation Set-up Swenson Walker Crystallizer Absorption in Packed Bed Liquid Liquid Extraction in Packed Bed Sieve Plate Distillation Column Spray Extraction Column Absorption in Sieve Plate Column Fluidized Bed Dryer Vapor in Air Diffusion Apparatus Forced Draft Tray Dryer Fixed Bed Absorption with Regeneration Solid in Air Diffusion Apparatus Natural Draft Tray Dryer York Scheibel Extraction Column Liquid in Liquid Diffusion Apparatus Rotary Dryer Batch Crystallizer Mass Transfer with/without Chemical Reaction (Solid-Liquid System) Solid-Liquid Extraction (Bonnoto type) Vapor-Liquid Equilibrium Set-Up Wetted Wall Column, Humidification & De-Humidification Study Experimental Water Cooling Tower Ion Exchange Mass Transfer Lab Equipment, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

