

Product Name :
Automatic Microtome Cryostat**Product Code :**
NMCBS13006**Description :**

Automatic Microtome Cryostat Our company caters to a wide market for Automatic Microtome Cryostat. These products are provided with solid state temperature controller -cum-digital indicator and is maintained with in $\pm 2^{\circ}\text{C}$ of set temperature. Automatic Microtome Cryostat was incorporated automatic defrosting system and delay timer. These are supplied completely with rust proof in build micro tome forced air refrigeration unit with transparent door cover, four demountable specimen holders, one bottles of low temperature oil and one micro tome knife 120 mm long with back and handle in wooden box, cord and plug. Automatic Microtome Cryostat Suitable to work on 220V, single phase, 50 Hz, AC supply. Control panel includes the following: • Digital Temperature Indicator cum controller • Defrost Time • Main switch • Delay Timer and cooling Indicator • Voltmeter

Optional: • Voltage Stabilizer (3 KVA) • Block Holder • Low Temperature oil Automatic Microtome Cryostat, Microtome Instrument Exporters, Microtome Instrument Manufacturers