



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Suction Machine

Product Code :
NMCBS16004



Description :

Suction Machine · One autoclavable polycarbonate jar with overflow valve system · One antibacterial and hydrophobic filter · One complete set of silicone tubes · One conical connector for cannula ct · One Power Cord
Technical specifications: · Motor: Oil less and maintenance free piston pump · Power Feeding: 220 50Hz other configurations on request · Maximum Vacuum (adjustable): -0.80bar -80kPa -600mmHg · Max free air Flow rate: 30l/min · Years of warranty: 1 · Shipping Carton: 4 Suction Machine, Suction Unit Exporters, Suction Unit Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,