

**Product Name :**  
Standard Consumables Suction Unit**Product Code :**  
NMCBS16005**Description :**

Standard Consumables Suction Unit Two autoclavable polycarbonate jars of 2000ml with overflow valve system  
· One antibacterial and hydrophobic filter · One complete set of silicone / PVC tubes · One conical connector for cannula · One Power Cord Options: Available with Oil Free Piston Pump Technical specifications: · Motor: Rotary Pump · Power Feeding: 220V 50Hz · Other configurations on request Maximum Vacuum (Adjustable): 700mmHg · Max free air flow rate: 50l/min · Duty Cycle: Non - stop operating · Years of warranty: 2 · Shipping Carton: 1 Standard Consumables Suction Unit, Suction Unit Exporters, Suction Unit Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,