



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Portable Suction Machine

**Product Code :**  
NMCBS16006



**Description :**

Portable Suction Machine Two autoclavable polycarbonate jars with overflow valve system · One antibacterial and hydrophobic filter · One complete set of silicone / PVC tubes · One conical connector for cannula · One Power Cord Options: Also available with Rotary Oil Dip Pump Technical specifications: · Motor: Oil less and maintenance free piston pump · Power Feeding: 220 50Hz other configurations on request · Maximum Vacuum (adjustable): -700mmHg+/-10mmHg · Max free air flow rate: 30l/min · Years of warranty: 1 · Shipping Carton: 1 Portable Suction Machine, Suction Unit Exporters, Suction Unit Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)