

**Product Name :**  
Digital Ultrasonic Cleaner

**Product Code :**  
NMCBS17002



**Description :**

Digital Ultrasonic Cleaner and Ultrasonic Cleaner for Hospital Lab and Medical Lab from India

**Technical Specification :**

**Digital Ultrasonic Cleaner Specifications:**

Tank Dim.: 45cm x 15cm x 15cm Liter

Cap.: 10 Ltrs U/S

Power: 200 Watts U/S

Freq.: 25 khz

Input: 180 Vac - 270 Vac

S.M.P.S. ensures smooth working of cleaners within the voltage variation of 180-270 vac.

Half the weight & size of the conventional machines

Heater ensures high efficiency in cleaning Cleaners are equipped with Imported Transducers for smooth and long life working.

Models can be made as per your specific requirements as per Industries Applicable, Hospitals, Nursing Homes, Educational Institutes, Electronic Industry, Jewelry Industry, Hotel Industry, Pharmaceutical Industry &

---

Automobile Industry **Digital Ultrasonic Cleaner, Ultrasonic Cleaner Exporters, Ultrasonic Cleaner Manufacturers**

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,