

Product Name :
Blood Bank Refrigerator**Product Code :**
NMCBS20003**Description :**

Blood Bank Refrigerator We have gained expertise in providing our customers with Blood Bank Refrigerator. Mounted on a sturdy angle frame, the inner chamber of our refrigerator is made of stainless steel and outer chamber of M.S. duly powder coated paint. In these products, temperature is maintained with 4 degree to 6 degree calcium for storing blood bottles or bags. As a safety device, our refrigerators are fitted with a power supply failure alarm. In addition to this, these products can be customized as per the specifications laid down by the customers. Other information: Provided with drawers for holding the bags to work 220-230V and 50 to 60HZ capacity Available in different capacities; 50 blood bags, 110 blood bags, 200 blood bags, 250 blood bags and 360 blood bags Blood Bank Refrigerator, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,