



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Plasma Extractor

Product Code :
NMCBS20011



Description :

Plasma Extractor Plasma Extractor is designed to Extract Blood component from centrifuged bags . The spring loaded front panel applies pressure on the collected bag causing the liquid to contain in transfer bag . Its sturdy construction makes it durable and reliable, easy to sue and portable . The unit is available in both versions automatic and manual Technical Specification: Sensor - Infra Red Clamping Control Manual Motor Activated Pressure Plate Acrylic Acrylic Alarm - Audio Visual Plasma Extractor, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)