

**Product Name :**  
Sphygmomanometer Aneroid Round Shaped - Mobile**Product Code :**  
NMCBS22002**Description :**

Sphygmomanometer Aneroid Round Shaped - Mobile ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Round Shaped - Mobile, OPD Instrument Exporters, OPD Instrument Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,