

Product Name :
Dual Bottle Suction Machine**Product Code :**
NMCBS24004**Description :**

Dual Bottle Suction Machine Each suction method independent of other · Comes with separate vacuum gauge support with reading from 0 to -760 mmHg · Electric unit option coming fitted with suitable motor support that provides for continuous flow rate · Unit comes with provision of foot operation in case of power failure · Bottle coming fitted with automatic float cut-out so as to protect pump from accidentally sucking in fluids · Front panel comes with support of pilot indicator lamp, ON/OFF switch, vacuum gauge and vacuum control knob Dual Bottle Suction Machine, Suction Apparatus Exporters, Suction Apparatus Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com