



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Horizontal Cylindrical Sterilizer

**Product Code :**  
NMCBS25002



**Description :**

Horizontal Cylindrical Sterilizer Steam Working Pressure 1.26 kgf/cm<sup>2</sup> (2.2 kgf/cm<sup>2</sup> in case High speed sterilizer  
Steam Working Temperature 121oC Operating Voltage 240, single phase, AC supply, 50 Hz · Sterilization  
Period: 45 to 50 minutes · Radiant locking system, cold water drainage system · Steam Exhaust: 5 to 7 provided  
with timer from 1 to 60 Min with double safety valve, steam releasing valve, pressure gauge for the pressure  
inside the chamber, and one pressure gauge to show the pressure of outer chamber both the chamber made of  
stainless steel lid of the unit made of heavy duty gun metal piece for extra safety, gasket, water level indicator to  
see the water inside, heavy duty power plug with socket provided with the unit, accessories: one spare element,  
and two gasket , our common size 400x1200 mm Horizontal Cylindrical Sterilizer, Medical Sterilizers Exporters,  
Medical Sterilizers Manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,