

**Product Name :**  
Digi Doppler**Product Code :**  
NMCBS26004**Description :**

Digi Doppler Table top · Rust proof · ABS Plastic molded · Rechargeable battery · Mains operated as well · Bright LED Display to Indicate the Beats per minimum (BPM) and the digital display of volume level · FHR with HI/Low limit setting · Audio Visual alarm Digi Doppler, Infant Care Products Exporters, Infant Care Products Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,