

Product Name :
Infant Baby Incubator**Product Code :**
NMCBS26005**Description :**

Infant Baby Incubator The temp controlled unit is operated by advanced micro computer technology · The three digital display indication provide the Set Temperature, Air Temperature and Skin Temperature · The safety features are self test and alarm system. Which will detect the different failures · The removable heater and fan unit make it easy to clean and service · In addition, the Skin/Air temperature display is larger then the other set temperature display · Full visibility & accessibility of the baby is provided by the designed of acrylic hood · A full complement of alarm provides safety convenience, with both audible and visual indicator Electrical Specification: Power 200V/50Hz Power Consumption 300Watts Chassis Leakage >500 Amp Current Fuse 5Amp Infant Baby Incubator, Infant Care Products Exporters, Infant Care Products Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,