



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Double Surface Phototherapy Unit

Product Code :
NMCBS26007



Description :

Double Surface Phototherapy Unit · Baby Bassinet with fold down panels · Head down position, on both sides · Phototherapy Unit fitted at a distance of 300 mm from baby bed · 4 Blue + 2 white tubes in Overhead Phototherapy Unit Choice of two kinds of tubes for use with · 4 nos. of blue light tubes with higher intensity of light in therapeutic range and 2 nos of white light tubes · 6 nos. of white light tubes with the facility to use 4 or 6 tubes at a time · Compact & mobile



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,