

**Product Name :**  
Cardiac Monitor SSM Cardiotrace

**Product Code :**  
NMCBS28001



#### Description :

Cardiac Monitor SSM Cardiotrace Tender Suppliers Africa and Cardiac Monitor SSM Cardiotrace Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Cardiac Monitor SSM Cardiotrace Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Cardiac Monitor SSM Cardiotrace China, Cardiac Monitor SSM Cardiotrace Kenya.

Cardiac Monitor SSM Cardiotrace Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Cardiac Monitor SSM Cardiotrace Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the

---

Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point), Suriname (Paramaribo)

#### **Technical Specification :**

Cardiac Monitor SSM Cardiotrace Owing to our experience, we have been successful in catering to the requirements of clients by offering quality Cardiac Monitor SSM Cardiotrace. These products are widely acclaimed by the customers for their excellent performance, easy installation, low maintenance and longer service life. In addition to this, we provide customization as per the specifications laid down by the customers.

Features : Cardiotrace is a portable CRT based Two Trace ECG monitor It is self contained unit for continuous ECG monitoring and display of heart rate Cardiotrace uses press button switch on the front panel ensuring trouble free operation Cardiotrace display 4 second of ECG on top channel and another 4 secs on the bottom channel Cardiotrace totally isolates the patient from the mains using special optical isolation technique Specially designed for Indian conditions (heart, humidity, voltage fluctuation) rugged, continuously operable & user friendly Second channel can be frozen of the detailed analysis. (It freezes automatically in the event of alarm violation) Bright 7 segment LED of 14 x 8 mm display for heart Rate which can be easily read from a distance also Adjustable alarm limits 0-250, (hi and low limit) stop and run facility for detailed analysis Cal for 1 mv pulse calibration Lead selector for selecting lead I, lead II lead III Defib protection : QRS indication The speed of the ECG waveform can be selected 25mm per second or 50mm per second 8 second ECG display Central Monitoring System for SSM Cardiotrace : Can be used in ICU, ICCU, AMC & Post operative ward Eight Bedside monitors can be interfaced User can change the 16 ECG colors & back ground colors 24 Hours heart rate trend User friendly through on screen help menu 30 sec real time ECG for the detailed analysis Individual HR alarm limit setting for each & every bed Individual speed control for bedside ECG waveform Date, Time, Patient Name, Age, Alarm Low Limit, Alarm Hi Limit Display Cardiac Monitor SSM Cardiotrace, medical monitoring system exporters, medical monitoring system equipments manufacturers

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,