

Product Name :
Double Wall Autoclave**Product Code :**
NMALV10003**Description :**

Double Wall Autoclave

Technical Specification :

Made of stainless steel with the inner chamber of a thick pressure gauge.
The outer surface has a mirror polish.
The lid of this autoclave is thick along with a safety valve, pressure release valve, and drain valve.
Spring loaded valve to adjust the pressure.
Works on 220V AC.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com