

**Product Name :**  
Stereo Acuity Test Kit

**Product Code :**  
NMCBS36002



**Description :**

Stereo Acuity Test Kit ½ Designed to rapidly test for amblyopia and strabismus using gross thru fine Stereopsis (4800 to 20 seconds of arc) ½ Graded circle test from 400 seconds now down to 20 seconds ½ No monocular clues with new technology ½ New improved booklet 1005 ½ Random Dot Stereopsis Test with LEA Symbols: ½ Designed to rapidly test for amblyopia and strabismus in early and non-readers and non-verbal children & adults ½ Expanded Random Dot LEA Symbols Test (500, 250, 125, 63 seconds of arc) ½ Graded circle test now down to 12.5 seconds with no monocular clues ½ New improved booklet 1010 ½ Butterfly Stereopsis Test with LEA Symbols: ½ Designed to rapidly test for amblyopia and strabismus using gross through fine Stereopsis (2500 to 20 seconds of arc) ½ Graded circle test from 400 seconds now down to 20 seconds with no monocular clues ½ New improved booklet Stereo Acuity Test Kit, Medical Keratometer exporters, Medical Keratometer manufacturers