

**Product Name :**  
Fully Automatic Autoclave**Product Code :**  
NMALV10004**Description :**

Fully Automatic Autoclave

**Technical Specification :**

Space between boiler &amp; Outer Shell works as an Air Insulation.

All system hydraulically tested upto 40 p.s.i as a safety measure Fitted with a (Programmable Logic Controller)

LCD Display with respect to time and temperature programmable by user enables automatic purging of stale air, sterilizing hold time for sterilizing period and cycle end with automatic exhaust of pressure.

Provided with Class 'A' Platinum (Pt-100) sensor, enabling precise and accurate monitoring &amp; Solenoid Valve for automatic purging and exhaust (user settable).

At the end of the cycle an audio alarm/ buzzer is raised in order to alert the operator.

