

**Product Name :**  
Vital Sign Monitors**Product Code :**  
NMCBS38005**Description :**

Vital Sign Monitors Vital Sign Monitors are 7"Color TFT Screen, waveform up to 8 channelsLight and portable with build-in rechargeable battery72 hour storage and review of trend grams and tables · Vital Sign Monitors are 40 second review of olographic wavefo RMS · Vital Sign Monitors are All-round monitor of adult, pediatric and neonatal electrosurgical unit and defibrillation protected · Vital Sign Monitors are Standard parameters: ECG, NIBP, SpO2, Respiration, Temperature, Pulse rate Vital Sign Monitors, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,