



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
6 Channel ECG Machine

Product Code :
NMCBS40002



Description :

6 Channel ECG Machine ½ Heart rate ½ Heart rhythm ½ For conduction abnormalities ½ Prior heart attack ½ Coronary artery disease ½ Abnormally thickening of heart muscle Features: ½ 6 Channel ECG Machine are 7" TFT, touch screen, led backlight ½ 6 Channel ECG Machine are Digital isolation technology and signal processing solution, digital filter ½ 6 Channel ECG Machine are Patient management, name and age and id. ½ 6 Channel ECG Machine are Detail analysis report. ½ 6 Channel ECG Machine are Multi-language options: Germany, France, Italian, Turkey, etc ½ 6 Channel ECG Machine are Built-in lithium rechargeable battery, wide thermal printing system ½ 6 Channel ECG Machine are USB and LAN socket. ½ 6 Channel ECG Machine are Memory: Built-in memory or mini SD Memory card. Store more than 1000 pieces of archive ½ 6 Channel ECG Machine are Auto-measurement, auto-analysis and auto-interpretation 6 Channel ECG Machine, ECG Machines exporters, ECG Machines manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,