

Product Name :
LP10 Ventilator**Product Code :**
NMCBS42009**Description :**

LP10 Ventilator
• Pressure control knob for pressure-limited, time-cycled ventilation
• Decreased internal resistance reduces breathing effort in SIMV mode
• Continuous, low-maintenance operation with minimal supervision results in low cost of ownership
• Features and portability provide flexibility to meet a broad range of ventilation and patient requirements
• Internal and external battery operation provides convenience and security
Indications For Usage:
• A portable ventilator ideal for in hospital transport, sub acute facilities and long term, in-home care

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com