

**Product Name :**  
Hand Held Pulse Oximeter**Product Code :**  
NMCBS43004**Description :**

Hand Held Pulse Oximeter Hand Held Pulse Oximeter are Digital Technology Hand Held Pulse Oximeter are Alarm with freely set limit Hand Held Pulse Oximeter are Accurate measurement during low perfusion Hand Held Pulse Oximeter are Anti-Movement, Stable Performance Hand Held Pulse Oximeter are Built-in rechargeable battery and seat-type charger, Power-Saved Function Hand Held Pulse Oximeter are Built-in memory card for data memory, connection with computer for downloading data Hand Held Pulse Oximeter are Bright, Easy-to-Read LED Numerals Hand Held Pulse Oximeter are Light weight, Portable 2 parameters: SPO2 Pulse Rate Hand Held Pulse Oximeter, Pulse Oximeter exporters, Pulse Oximeter manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)