

**Product Name :**  
Table Top Pulse Oximeter With NIBP

**Product Code :**  
NMCBS43005



**Description :**

Table Top Pulse Oximeter With NIBP Table Top Pulse Oximeter with NIBP are Displays SpO2, Pulse Rate, and Plethysmograph Table Top Pulse Oximeter with NIBP are Digital SpO2 Technology to reduce alarms and improve performance during motion and low perfusion Table Top Pulse Oximeter with NIBP are Graphical & Tabular Trend Memory Table Top Pulse Oximeter with NIBP are Excellent Low Perfusion Performance Table Top Pulse Oximeter with NIBP are Audible and Visual Alarms Table Top Pulse Oximeter with NIBP are Over 4 hours battery backup Table Top Pulse Oximeter with NIBP are Automatic Battery Charging Table Top Pulse Oximeter with NIBP are SpO2 dependant pulse tone Table Top Pulse Oximeter With NIBP, Pulse Oximeter exporters, Pulse Oximeter manufacturers

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)