



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
100MA Mobile X Ray Machines

Product Code :
NMCBS47004



Description :

100MA Mobile X Ray Machines Counter-balanced. tube arm for increased reach better positioning Integrated cassette box for easier operation and better cassette handling Elegant digital display and extremely easy to use controls MA meter & voltmeter for line and exposure monitoring Automatic calculation of mass with digital display Digital display of KVP 6 Led indicator for line, 0/L, X Ray, Stand by, Bucky & fuse Precision solid state electronic time.



Medical Hospital Equipments Manufacturers

Website: SITE