

Product Name :
Digital Hemoglobin Meter

Product Code :
NMCBS58002



Description :

Digital Hemoglobin Meter Instrument Digital Hemoglobin Meter Display Digital 3½ Digit, Bright LED Display Mode % T, O.D. and 3 concentration Modes, C1,C2,C3 Zero Adjustment Automatic (Without Knob rotation) Accuracy + 0.02 O.D. Stability + 0.02 O.D. Per Hour Light Source 6.8 Volt, 0.3 Amp Tungsten Lamp Amplifier Micro controller Based Memory In built memory to retain last 1½ factor for concentration even after power loss. Sample Volume 1 ml. minimum Power 230 V + 10 % AC, 50 Hz. Accessories Set of Test Tube- 5 No. Dust Cover 1 No. Operation Manual Salient Features: 1½ Fully Automatic 1½ Digital Read Out 1½ Highly Accurate 1½ Long Life Digital Hemoglobin Meter, Hemoglobin Meter exporters, Hemoglobin Meter manufacturers