

Product Name :
Haemoglobinometer**Product Code :**
NMCBS58003**Description :**

Haemoglobinometer It is black in colour. The back of the holder is fitted with milky colour screen made of acrylic. In front of the holder there are three windows. The two side windows are fitted with colour comparison. The size of all the three windows are Height 3 cms. (Minimum) Width 0.5 cms. (Minimum) ½ Square H.B. Tube: -2 Nos. ½ Square H.B. Tube is graduated on both the side for ½ the measurement of haemoglobin in respect of ½ percentage and gram. ½ H.B. Pipette -1 Nos. ½ 20 microlite H.B. Pipette with latex tube and mouth piece ½ Amber coloured glass bottle -1 Nos. ½ Glass Stirrer -1 Nos. ½ Dropper with teat -1 Nos. ½ Cleaning Brush -1 Nos. ½ Spare latex tubing one meter -1 Nos. ½ All the above are fitted in thermoform tray. Further packed in plastic box. Haemoglobinometer, Hemoglobin Meter exporters, Hemoglobin Meter manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,