

Product Name :
Portable Semi Automatic**Product Code :**
NMALV10013**Description :**

Portable Semi Automatic

Technical Specification :

Top cover clamped to the lower tank by wing nut closing system

The unit has a spring-loaded safety valve for pressure adjustment and a dead weight safety valve. Supplied with electric fitting, cord, and plug.

Top cover seals the autoclave with neoprene rubber gasket

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,