

**Product Name :**  
Portable Semi Automatic**Product Code :**  
NMALV10013**Description :**

Portable Semi Automatic

**Technical Specification :**

Top cover clamped to the lower tank by wing nut closing system

The unit has a spring-loaded safety valve for pressure adjustment and a dead weight safety valve. Supplied with electric fitting, cord, and plug.

Top cover seals the autoclave with neoprene rubber gasket

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)