



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
DC Motor with Techo Generator

Product Code :
NMCBS88016



Description :

DC Motor with Techo Generator DC Motor with Tacho generator. Being a quality-oriented firm, we assure for the quality and flawlessness of this motor. This motor is extensively used in power plants, chemical plants and other industries for various purposes. Moreover, our offered range can be availed in various technical specifications, as per the requirements of our esteemed clients. Features: DC Motor with Techo Generator are Low output rotation speed DC Motor with Techo Generator are low noise block design DC Motor with Techo Generator are Allow multi-stage combination DC Motor with Techo Generator, Ac Motor exporters, Ac Motor manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,