

**Product Name :**  
Infant Radiant Warmer, Open Care System

**Product Code :**  
NMBCE10003



**Description :**

Infant Radiant Warmer (Infant Resuscitation Trolley), Open Care System on wheels. Consists of Ceramic Infra Red Heater 650 Watts. Metallic overhead heater box. The construction of the stainless steel reflector and guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards. Halogen Examination Lamp (12 Volts, 50 Watts). To provide accurate assessment of Infant colour and sufficient illumination during procedures.

Essential Attachments: S.S. IV Stand. Infant Care Trolley Standard. Clear acrylic collapsible side panels to assess an infant from all sides. Stainless steel baby tray with transparent slot to incorporate the facility of under surface photo-therapy. Acrylic tray and cushion mattress with head up/down facility. Acrylic X-ray tray with X-ray guide for positioning the cassette. Two Lower shelves. Mounted on wheel 2inch, The design aspects also allow easy of dismantling baby care area during cleaning or fumigation.

Feature: Control Modes: Skin / Air / Manual

Range: Skin control: 34 to 38 C Air control: Room temp. to 39 C

Resolution: 0.1 C

Accuracy:  $\pm 0.2$  C

Display: 1in bright skin & air temp.(7 segment LED)

Status: On LCD

Control Mode & Set temp

Timer: Timer / APGAR Timer / Baby Attend

Manual mode: Time Programmable

Heater output settable

Over temp. cut-off

---

Alarm: Auto audio alarm mute, Mute time programmable

Safety alarms: ON LCD

Heater Failure

Power Control Failure

Probe & Power failure alarm Infant Radiant Warmer, Open Care System, Baby Care Equipments Exporters,  
Baby Care Equipments Manufacturers

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)