

**Product Name :**  
Soil Testing Kit**Product Code :**  
NMCBS104003**Description :**

Soil Testing Kit pH Range 0 to 14pH Temp. Comp. 0oC (Manual) Resolution 0.01pH Accuracy +0.01 pH  
O.R.P./mV Range N.A. 0 to + 1999 mV Resolution 0.01pH N.A. 1 mV Accuracy N.A. + 1 mV + 1 digit  
Conductivity Range 1 5 Cell Constant 0.4 - 1.5 Adj Resolution 0.1 Micro mhos Accuracy + 0.5% FS Temperature  
Range 0o - 100oC Resolution 0.1oC Accuracy +0.1oC + digit T.D.S. Range 1 5 Cell Constant 0.4 - 1.5 Adj  
Resolution 0.1ppm Accuracy +5% FSD D.O. Range 0-20 ppm N.A. Temp. Comp. 0.1 ppm N.A. Resolution +0.3  
ppm N.A. Accuracy 5o to 50o (Manual) N.A. Turbidity Range 0-1000 NTU/JTU N.A. Resolution 1 NTU/JTU N.A.  
Soil Testing Kit, Digital Water & Soil Analysis Kits exporters, Digital Water & Soil Analysis Kits manufacturers

**Technical Specification :**

Automation Grade: Manual

Usage: To determine on the spot soil pH value. This is a compact kit easy to carry on site.

Color: Other Color

Soil Elements: Nitrogen, Potash, Phosphorus

No of Equipment: 12

Usage: Soil Test

