

Product Name :
Naugramedical Fowler Bed Manual**Product Code :**
NMCBS106002**Description :****Naugramedical Fowler Bed Manual India China****Technical Specification :**

Naugramedical Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism

Size: 218 940 540 mm.

Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward.

Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel.

Epoxy coated mild steel frame work and 4 section perforated top.

Circulated Central braking system locks / unlocks all swivel castor simultaneously.

Provision for I.V. Rod. Fowler Bed Manual, Hospital Bed exporters, Hospital Bed manufacturers

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]