

Product Name :
I.C.U. Bed, Electric, Six Function with Sitting Position

Product Code :
NMCBS106003



Description :

I.C.U. Bed, Electric, Six Function with Sitting Position & Comfy ICU Bed (6 function) with electrically operated sitting position 0-70°, back rest tilting 0 - 75°, knee rest tilting 0-35°, trendelenburg tilting 0 - 15°, Reverse trendelenburg tilting 0 - 15°, & Size: 2155 x 950 x 570 / 780 mm (Adjustable Height). All functions controlled with Electric Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed. I.C.U. Bed, Electric, Six Function with Sitting Position, Hospital Bed exporters, Hospital Bed manufacturers