

**Product Name :**  
ICU Bed Electric 7 Function

**Product Code :**  
NMCBS106005



**Description :**

**ICU Bed Electric 7 Function India China**

**Technical Specification :**

I.C.U. Bed, Electric, 7 Function

ICU Bed (7 function) with electrically operated back rest tilting 0-75

knee rest tilting 0-45

trendelenburg tilting 0-27

Reverse trendelenburg tilting 0-27

mattress base tilting to the left up to 3

tilting to the right up to 35

Size: L 2170mm

W900mm

H 520-910mm(Adjustable Height).

All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request).

---

Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward.

Easy to operate built-in Control Panel on both sides of guard rails.

Separate Movable Nurse Control Panel can lock the buttons of Guard Rail Control panels and Foot Step control panel.

High quality beautiful ABS engineering plastic head panel and foot panel equipped with bumper protection.

Epoxy coated mild steel frame work and 4 section perforated top.

Equipped with battery backup.

Noiseless castors with Central braking system.

Provision for I.V. Rod on both sides of the bed. I.C.U. Bed, Electric, 7 Function, Hospital Bed exporters, Hospital Bed manufacturers



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,